

OVER THE EDGE PATIENT PACKAGING SKILLS EVALUATION

Guidelines:

Rescuers packages patient appropriately based on conditions present and scenario given.
 Rescuers packages patient with equipment provided in separate caches.
 Rescuers to demonstrate proficient knowledge of packaging technique and ability to adapt.

Required Equipment:

- | | |
|--|---|
| C-spine head/neck | C-spine Spiderstrap Method |
| Head rolls/blocks | Spiderstraps |
| Tape | C-Spine Webbing Method (1" webbing) |
| Head pad (optional) | 15' Webbing—Chest harness |
| Flight Harness—chest and pelvic w/ carabiner | 15' Webbing—Pelvic harness |
| Stokes litter | 20' Webbing—Shoelace to backboard |
| One 25' webbing—external shoelace | Backboard |
| or two 15' webbing—external shoelace | Two 5' webbing—secures opposing corners |

Scenario: Patient with c-spine injuries under flyable conditions. Demonstrate the following patient packaging skills.

Observed Skills:

	Yes	No	Notes
Manual C-spine (verbalized)			
Simulate marrying head			
Log roll to backboard.....			
Attach Flight Harness (for flyable conditions)			
Proper chest placement and attachment.....			
Proper pelvic placement and attachment			
Attached carabiner and chest strap to pelvis.....			
Adjusts harness as packaging conditions dictate ..			
C-Spine to Backboard—Spiderstrap method.....			
Proper shoulder placement and attachment			
Proper pelvic placement and attachment			
Proper foot attachment and tension			
Aware of knee placement (bonus if verbalized) ..			
C-Spine to Backboard—Webbing method			
Proper knots (2 round turns & 2 half hitches)			
Seat harness—secured properly w/tension			
Chest harness—secured properly w/tension			
Shoelace patient to backboard			
C-spine head and neck			
Apply head padding (optional)			
Apply head blocks			
Apply tape.....			
Backboard to Stokes			
Secure stokes head and foot short webbing			
Shoelace webbing—proper tension and knots.....			